



## GB MMLV reverse Transcriptase

For Research use only

Cat No: GB2430

Size : 5 ml ( 1,000,000 units)

### Description:

GB M-MLV Reverse Transcriptase RNase H Minus is a modified version of the Moloney Murine Leukemia Virus Reverse Transcriptase, designed to remove the RNase H activity. This modification makes it particularly useful for synthesizing long cDNA from long RNA templates (greater than 5kb).

### General Protocol

- 1. Template RNA Preparation:** Isolate RNA from cells or tissues and ensure it is free of DNA contamination.
- 2. Reaction Setup:**
  - RNA Template: Typically 1 µg of total RNA or 10-100 ng of mRNA.
  - Oligo(dT) Primer or Random Primers: Primers to initiate cDNA synthesis.
  - dNTPs: Deoxynucleotide triphosphates for DNA synthesis.
  - M-MLV Reverse Transcriptase RNase H Minus: The enzyme for cDNA synthesis.
  - Reaction Buffer: Provides the necessary ions and environment for enzyme activity.
- 3. Reverse Transcription:**
  - Annealing: Incubate the RNA-primer mix at 65°C for 5 minutes, then cool on ice.
  - Extension: Add the M-MLV Reverse Transcriptase RNase H Minus and incubate at 37-42°C for 50-60 minutes.
  - Termination: Heat at 70°C for 15 minutes to inactivate the enzyme.

### Advantages

- **High Yield:** Produces high yields of cDNA, making it ideal for downstream applications.
- **Robustness:** Effective with various types of RNA, including challenging templates.
- **Ease of Use:** Simplifies the process of reverse transcription with standardized protocols.

### Limitations

- **RNA Quality:** Requires high-quality RNA to produce reliable results.
- **Enzyme Sensitivity:** Sensitive to inhibitors commonly found in RNA preparations, such as ethanol or salts.

GB M-MLV Reverse Transcriptase RNase H Minus is a valuable tool for molecular biology research, enabling the conversion of RNA into cDNA for various applications.

